

Subject: Mirth Connect Security Vulnerabilities

Date: 2024-11-29 (last update)

Background

NextGen Healthcare Mirth Connect is an open-source integration engine used primarily in healthcare IT for exchanging healthcare data between various systems. It enables interoperability between different healthcare applications, and allows secure and efficient transfer of data through standardized protocols and formats such as HL7, DICOM, and FHIR.

IHTeam identified a remote command execution on **Mirth Connect <= 4.3.0** that can be exploited from an unauthenticated perspective. The only condition in which this Common Vulnerabilities and Exposures (CVE) can be exploited is when **JRE <= 8.0** is in use. The vulnerability reportedly came as a side-effect of the company trying to fix a previous critical-severity flaw, tracked as [CVE-2023-37679](#). This vulnerability, carrying a severity score of 9.8, was also described as a pre-auth remote code execution.

Later, researchers from Horizon3.ai determined the patch to be incomplete and published a gadget chain which bypassed the deny list that the original had implemented. This second vulnerability was assigned [CVE-2023-43208](#) and was patched in **Mirth Connect v4.4.1**.

Response

Affected Vyairé Products

The Vyairé SentryConnect solution, which is utilized to connect SentrySuite™ solutions to hospital information systems, utilizes [Mirth Connect](#) for integration purposes. **Mirth Connect before v4.4.1** is potentially vulnerable to [CVE-2023-37679](#) and [CVE-2023-43208](#).

How Vyairé Is Responding

Vyairé has implemented the below changes around the disclosure of the Mirth Connect vulnerabilities from **SentryConnect Gateway version 6.3.1.7** onwards.

- Updated to **openJDK v17.0.11+9-LTS**
- Mirth server listening on **localhost only**
- **Removed insecure cipher suites:** TLS_RSA_WITH_AES_128_GCM_SHA256 and TLS_RSA_WITH_AES_256_GCM_SHA384
- Upgraded to **MirthConnectAdminLauncher v1.4.1**
- Upgraded to **MirthConnect 4.50**
- Upgraded MirthConnect's jetty 9.4.49 to **jetty 9.4.53**

Product Security Bulletin

Important Information - Please Read and Keep



Mitigations & Compensating Controls

Vyairé recommends applying the following mitigations and compensating controls:

- Customers using versions before **SentryConnect Gateway version 6.3.1.7** must reach out to **support.connect.eu@vyairé.com** to initiate the update process and also allow Vyairé support to make additional configurational changes.

For product or site-specific concerns, contact your Vyairé service representative.

For more information on the Vyairé proactive approach to product security and vulnerability management, contact us at productsecurity@vyairé.com or visit www.vyairé.com/product-security.