

Save time, maximize resources and improve workflow with Avea™ disposable expiratory filter/water traps

Disposable option

Step 1: Remove

Simply remove the Avea™ disposable expiratory filter/water trap from your AVEA ventilator and dispose of it.



Step 2: Replace

Replace the used filter with a fresh Avea disposable expiratory filter/water trap and you're ready for the next patient!

Reusable option

The single-patient-use Avea disposable expiratory filter/water trap incorporates a simple 2-step process that saves you time, maximizes your resources, improves workflow and reduces the risk of process variability.

Reusable filters and water traps may cost less up front but they consume hospital time, energy and resources across a multi-step tangle of personnel, departments, processes and protocols.

The multi-step workflow for processing of reusable filter and water trap components:



Comparison of benefits



Avea disposable expiratory filter/water trap

	Avea disposable expiratory filter/water trap	Reusable filter/water trap
2-step replacement process improves workflow efficiency	•	
Single-patient-use reduces risk of process variability	•	
No component assembly required	•	
Easier condensation management	•	
HEPA filtration element filters the expired gas to 0.3 microns	•	•
99.999% Viral & Bacterial Filtration Efficacy (VFE/BFE)	•	
60% lower flow resistance than the reusable configuration	•	
Accommodates neonate through adult patients	•	•

Avea disposable expiratory filter/water trap

Part number: 11790

Quantity: Case of 12

For more information, contact your local sales representative.

Additionally, please feel free to call Order Administration at **800.231.2466**.

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA

vyaire.com

REG ONLY For U.S. distribution only.

©2018 Vyaire. Vyaire, the Vyaire Logo and Avea are trademarks of Vyaire Medical, Inc. | CF/899/18/0001 (0318)

vyaireTM
M E D I C A L